

FILLOGY
FILLING EVOLUTION

FILLING KIT FOR COSMETS
AND DRUGSTORE ARTICLES





FILLOGY LIQUID XS, Art. 301043

Small amounts of liquid products, even chemical or abrasive products, e.g.

- perfumes
- facial oils
- dyes, flavours
- sterile products

Product feeding: via hose \varnothing 1,6 mm/3,2 mm/6,4 mm/8,0 mm
 Filling speed: 175 - 2,600 ml/min (depends on the hose- \varnothing)
 Filling amount: 2 ml - 1,000 ml
 Filling accuracy: +/- 1 %
 Filling temperature: 2 - 95 °C / 35,6 - 203 °F



FILLOGY FLOWABLE S, Art. 301019

Flowable products, also with microparticles, e.g.

- fingernail UV gel
- polishing gel
- lotions
- hair colours

Product feeding: via funnel of 5 l
 Filling speed: 800 ml/min
 Filling amount: 2 g - 100 g
 Filling accuracy: +/- 1 %
 Filling temperature: 2 - 95 °C / 35,6 - 203 °F



FILLOGY CREAMY S, Art. 301058

Creamy products, also small amounts, e.g.

- cremes
- balms
- body lotions
- mascara
- make up

Product feeding: via container of 3 l
 Filling speed: up to 5 l/min
 Filling amount: 2 g - 100 g
 Filling accuracy: +/- 1 %
 Filling temperature: 2 - 95 °C / 35,6 - 203 °F



AUTOMATABLE WITH SMALL TURNTABLE FILLOGY AUTO S, Art. 308174

For semi-automatic filling of container sizes from 2 ml to 500 ml
especially suitable for containers with small opening diameter

Rotation speed: 0,4 - 4,5 U/min adjustable
 Diameter of container opening: $d \geq 8,0$ mm
 Storage space: 35 x 500 g glasses (15 empty + 20 full)
 45 x 250 g glasses (19 empty + 26 full)
 Diameter turntable: 635 mm
 Load-bearing capacity: 20,0 kg
 Footage: Width 635 mm x Depth 870 mm
 Total Height: 655 mm / 1.385 mm with rolling frame

